



Windows 10 ABC Client Laptop Setup for Remote Shadow RDP

Contents

Configuring the client laptop when it is switched on for the first time.....	3
Section 1 Configuring the client laptop when it is switched on for the first time	3
WARNING 1.1	3
Section 2 Configuring the user accounts	3
Section 3 Configuring the computer power settings	4
Client (Laptop) configurations	5
Section 1 Client IP Configurations.....	5
Section 2 Setting up Client Shadow RDP Configuration.....	10
Installing USB over Ethernet on the Client side	10
Setting up the New RDP Icon on the Client PC	15
Testing Shadow RDP	15
Testing Client Connection with Gating Kit	17

Configuring the client laptop when it is switched on for the first time

Section 1 Configuring the client laptop when it is switched on for the first time

This is the procedure to configure the Active Breathing Coordinator™ computer when it is switched on for the first time.

If it is necessary, you must configure the computer when it is switched on for the first time.



WARNING 1.1

Do not install software which is not approved by Elekta. The result can be system malfunctions. If you ignore this warning, you can cause clinical mistreatment.

- 1 Select the **Language, Country/Region, and keyboard layout**, click **Next**.
- 2 Connect the network to internet via Ethernet or WiFi, or click Use Express settings or if unable skip this step.
- 3 Accept the license agreement.
- 4 Select set up **Organization** and click **Domain join instead**, click **Next**.
- 5 Type a user name **ABC administrator**.
- 6 Type **elekta** as the password. Input security questions.
- 7 Accept **Privacy Settings**, skip **Customize experience**, and select Not now **Cortana help**.
- 8 Click **Next** for **Preferences** and check mark **Dell's Privacy Policy and Telemetry Notice**.
- 9 Select **Decline updatgrade** for **Windows 11**.

Section 2 Configuring the user accounts

This is the procedure to configure the user accounts on the Active Breathing Coordinator™ computer.

Configure a User account on the Active Breathing Coordinator™ computer Switch on the Active Breathing Coordinator™ computer.

- 1 Log on as the **ABC Administrator** user.
- 2 Configure the laptop with the ABC administrator and ABC user accounts (see table below).


User name	Type	Password
ABC administrator	Administrator	elekta
ABC	Standard user	abc

- 3 Open **Settings**, click **Accounts**, and select **Family & other users**.
- 4 Select **Add someone else to this PC**. A dialog box will open.
- 5 In the username text box, type **ABC**.
- 6 Type **abc** in the password text boxes. Input security questions.
- 7 Click **Next**, and verify the user is a **Standard user** by clicking the name and selecting **Change account type**.

Section 3 Configuring the computer power settings

This is the procedure to configure the power settings on the Active Breathing Coordinator™ laptop. The laptop power settings must be correctly configured.


- 1 Log on as the **ABC administrator**.
- 2 Open the **Control Panel**.
- 3 Make sure that **View By** is set to **Category**.
- 4 Select **Hardware and Sound**.
- 5 Select **Power Options**.
- 6 Click **Change plan settings** under **Balanced**.
- 7 Set **Turn off the Display** to **Never**.
- 8 Set **Put the computer to sleep** to **Never**.



- 9 Click **Save changes**.


Client (Laptop) configurations

Section 1 Client IP Configurations



Step 1.

- Left Click **Internet Icon** “the globe Icon” or the **WiFi** signal icon.



Step 2.


Active Breathing Coordinator™ R3.0
Installation and Service

Copyright ©2023 Aktina All rights reserved.

Document ID20230112


Page 5

- Left Click “**Network & Internet settings**”




Step 3.

- Left click “**Change adapter option**”




Step 4.

- Double left click “**Ethernet**”




Step 5.

- Double left click **“Internet Protocol Version 4 (TCP/IPv4)”**




Step 6.

- Switch Radio Button to **“Use the following IP address”**




Step 7.

- Fill in the IP Address: **“192.168.1.55”** Subnet Mask:**255.255.255.0**
- Click **Ok**.




Step 8.

- Click **“OK”** to close



Step 9.

- Close the window by hitting the X on the upper right.





Step 10.

- Close the window by hitting the X on the upper right


Section 2 Setting up Client Shadow RDP Configuration Installing USB over Ethernet on the Client side

Step 1.



- Press the Windows Key and start typing "Edge"
- Press enter or click "Open" on the right-hand side to open "Microsoft Edge".
- While connected to the Internet or downloading from another PC with internet access. Go to the following URL link: <http://www.usb-over-ethernet.com/downloads/usb-over-ethernet.zip>
- Hit Enter, once typed into a web browser, it should automatically download a .ZIP file.
- Select "Open file" once its download.

Step 2.




- Once open, Double click “Setup”
- Follow the prompts to install the software.
- Select “I accept the agreement”
- Then click “Next” or “Install” while leaving the settings to default.
- Until you get to “Finished”
- Click “No” if it asks to perform automatic new version.

Step 3.




- Once the application starts.
- Locate the “Elekta Response Gating System (Communication Modem Device)”
- Please make sure the Gating kit is connected to the Client’s Laptop.
- Right click Elekta Response Gating System and Select “Share”

Step 4.




- Finish up by Selecting “Settings” on the top.
- Select “Preferences...”
- Under “General” tab make sure the first 4 boxes are selected and the last one turned off.
- Select “OK” once you finished.

Step 5.



- Not with the application still open.
- Click on “Help” and select “Enter License Key...”
- This will open a window that will allow you to copy the Product or manually type it in. You will need a license key for the USB over Ethernet install. Email support@aktina.com for permanent product key (a 15 days trial period is in effect at installation).
- Once entered you can Select “OK”
- If done correctly you will get a Prompt asking to restart the application.
- We recommend doing a restart on the device and log into the account that will mainly be used.


Step 6.



- By bring up the application once more.

- You may check to make sure its activated if logged into another profile.
- Navigate to “Help” on the top and select “About”
- Close out by clicking “OK”

Step 7.



- Lastly you can disable the sharing of any other device to mitigate confusion.
- Right click on the device you wish not to share and select “Properties”
- Make sure you check the box for “Deny manual and automatic sharing of this device”
- You’ll see a red icon device to show its disabled.
- Click “OK” to close.

Setting up the New RDP Icon on the Client PC

Step 1

Log into the ABC on Client Laptop.

Username: ABC

Password: abc

Step 2

Open "Microsoft Edge"

Go to url link and type <https://www.aktina.com/products/elekta-abc-system/abc-downloads/software-and-downloads/>

to download the latest files required for this instruction set. There will be one file to download: "ABC-Dual-Monitor.zip"

Step 3

Open ABC-Dual-Monitor.zip

Open ABC-Dual-Monitor-Setup\ABC-Client folder and copy ABC.bat file to the laptop's desktop. There should be the ABC icon on the desktop.




Note: There will be two directories in the zip file: "ABC-host" and "ABC-client".
The "host" is the computer within the ABC cart (**Not Needed for Laptop configuration**)
The "client" is the laptop computer that is connected to the host via Remote Desktop.

Testing Shadow RDP

Step 1


- Connect laptop to ABC Host via ethernet cable.
- Turn on ABC Host machine and log on as **ABC administrator**, password **elekta**
- Double click on "**Init**" icon desktop to reboot system.
- Log on to the ABC Host machine as **abc**, password **abc**

Step 2



- Double click the file ABC.bat on the laptop. Note: A dialog box "Windows protected your PC" may show for the first time, click "More info" to run the ABC.bat file.
- Type ABC on the top box and type abc as the password and click "OK" or press enter

Step 3




- You're now able to control ABC Host while still allowing video output on the monitor
- If you get an error because of the session ID or shadow ID. Reboot the ABC host and log in as ABC, password abc using the Restart shortcut that was made on ABC profile desktop

Testing Client Connection with Gating Kit

Verifying Shadow RDP connection.

Note: Make sure Gating Kit is connected to the laptop (via USB cable) and turned on. Ensure ABC Control module is on and ABC enclosure and patient switch are connected.


Step 1 Click ABC.bat to start connection to ABC Host



Step 2

Open ABC 3.0 on host machine

Navigate to Help – Response Diagnostics



Verify the following items:



- Communication Link UP
- Enabled (Note: if not enabled, make sure Gating Kit is on and connected to the laptop)
- Relay Link Up
- Open