



ABC Remote RDP Patch

Addendum to ABC Installation

Table of Contents

Configuring the Existing ABC Host to use RDP-Shadow.....	2
Setting up the New RDP Icon on the Client PC	13
Installing USB over Ethernet on the Client Side	22
Installing USB over Ethernet Client on the Host	27

Configuring the Existing ABC Host to use RDP-Shadow

Step 1


Go to

<https://www.aktina.com/products/elekta-abc-system/abc-downloads/software-and-downloads/>

to download the latest files required for this instruction set. There will be one file to download: "ABC-Dual-Monitor.zip"

Step 2

Log into the ABC on-board PC.



Username: ABC Administrator

Password: elekta

Step 3

Use a USB stick to transfer the downloaded file to the computer ABC onboard computer.

Step 4

Extract the ZIP file to:


C:\users\public\desktop

There will be two directories in the zip file: "ABC-host" and "ABC-client".

The "host" is the computer within the ABC cart


The "client" is any computer that is connected to the host via Remote Desktop.

Step 5



- Press the Windows Key and start typing “Task Scheduler”
- Left click “Open” on the right-hand panel.


Step 6



- With Task Scheduler open
- On the right-hand panel select “Import Task”
- Navigate to C:\Users\Public\Public Desktop by using the left-hand panel.
- To locate “ABC Banner” file.
- Select the file and left click “Open”


- Verify that under Security Options you see ABC Administrator. If you do, click OK to finish the import and skip to Step 8

Step 7




- If the name is a string of random number and letters.
- Select "Change User"
- In the new windows type "abc administrator" and left click "Check Names"
- Windows will auto fill the rest of the information.
- Click OK once filled in.
- Check to make sure the User account now has "ABC Administrator"
- Left click "OK" to finish the import.

Step 8




- Locate the new task created "ABC Banner"
- With the task highlighted, On the right-hand panel left click "Run" to test the Banner.
- Left click "OK" to close the window.
- Close out and return to the desktop. Skip to Step 16
- If banner failed cancel the import, proceed with the following steps to manually configure the task.

Step 9




- While still in the Task Scheduler, on the right-hand panel. Select “Create basic task”
- Type the following into the fields - Name: ABC, Description: abc warning message. Click “Next”
- Select “When I log in” and click “Next”

Step 10




- Select “Start a program” and click “Next”
- Select “Browse” to open file explorer – Navigate to C:\Windows\System32\WindowsPowerShell\v1.0
- On the left-hand panel select “OS C:” Drive

Step 11




- Double click and navigate the following folders “Windows”, “System32”, “WindowsPowerShell”, and “v1.0”.

Step 12




- Select “powershell” and click “Open”
- Add the following in the “Add argument (optional)”
- Type it out or follow the next steps to be able to copy it from the failed ‘ABC Banner.xml!’
- -WindowStyle hidden -Command "& {[System.Reflection.Assembly]::LoadWithPartialName('System.Windows.Forms'); [System.Windows.Forms.MessageBox]::Show('!!!!!!Reminder!!!!!! Logon as ABC user and run INITABC after doing any maintenance .')}"

Step 13




- Locate the “ABC Banner” file that was downloaded and extracted in steps 2 – 4.
- Right click the file and select “edit”
- Scroll down to the end of the file to locate `<Arguments>`
- Select the following text that’s between `<Arguments>... </Arguments>` as shown in the example.
- Right click and copy the string.

Step 14




- Type or Paste the String in the “Add Argument (optional)” and click “Next”
- Verify and click “Finish” to create the task.

Step 15




- On the left-hand panel double click on “Task Scheduler Library”
- Locate and select “ABC Banner” to highlight it
- Then select “Run” to test the Banner
- Click okay to close the “ABC Banner”
- Close out of any open windows to return to the desktop

Step 16




- From the desktop, double left click “Shadow Firewall Rule”
- If prompted by Open File _Security Warning, click “Run”
- Left click “YES” on the first window prompt and then “OK” on the second one.

Step 17



- Press the Windows Key and start typing “Windows Defender Firewall”
- Press enter or click “Open” on the right-hand side




- On the left-hand side Double click “Advanced settings”

Step 18



- Double click “Inbound Rules” on the left hand-side
- This should show all current rules.

Step 19




- Search for “ABC” to make sure you see the rule created.
- Close out the window once you have confirmed.

Step 20

Name	Date modified	Type	Size
ABC Banner	8/9/2022 5:01 PM	XML Document	4 KB
InitABC	7/29/2022 9:26 AM	Shortcut	1 KB
SessionID	7/29/2022 7:50 AM	Windows Batch File	1 KB
Shadow Firewall Rule	8/4/2022 11:56 AM	Registration Entries	1 KB
ShadowRegistryEdit	8/9/2022 9:19 PM	Registration Entries	1 KB

- From the Desktop locate “ShadowRegistryEdit”
- Double click the file to execute
- If prompted by Open File _Security Warning, click “Run”
- Click “Yes” and “OK” just like step 16.


Step 21



- From the desktop press the “Windows Key” and Type “Registry Editor”


- Select “open” to launch the application
- Navigate the folders on the left-hand panel to reach the following:
Computer\HKEY_LOCAL_MACHINE\SOFTWARE\Policies\Microsoft\Windows NT\Terminal Services
- Double click “HKEY_LOCAL_MACHINE”, “SOFTWARE”, “Policies”, “Microsoft”, “Windows NT”, and Select “Terminal Service” to highlight it.
- On the right-hand panel you should see an entry for “Shadow”
- Once confirmed you can close this window and skip to Step 23
- If you don't see “Shadow” already made. Proceed to the next Step.

Step 22




- Right click anywhere in the right-hand panel.
- Using the mouse hover over “New” to access the subcategory
- Select and click on DWORD (32-bit) Value
- Type: Shadow and press enter to confirm.
- Double click on newly created “Shadow”
- In the “Value data” Type: 2
- Leave the rest alone and click OK
- Close Registry editor to return to the Desktop.

Step 23




- While back at the Desktop, double click on session ID to run the file.
- This is used to troubleshoot what session ID ABC user has been assigned.
- We always want this to stay at Session ID 1

Step 24




- Before sign into ABC user.
- You can now right click on "ABC Banner" and select "Delete"
- Do the same for "Shadow Firewall Rule" and "Shadow Registry Edit"

Step 25



- Now we can check to make sure all it's working correctly.
- While at the desktop Hit the Windows Key
- Left click on user icon on the left-hand panel.
- Select ABC user to sign in.
- Log into ABC
- Password = abc

Step 26




- Double click on INITABC icon
- Enter Admin credentials
- Password = elekta
- Wait and you should see that system will restart.
- Once the system has turned back on. Log into ABC

Setting up the New RDP Icon on the Client PC

Install setup:


- Connect the client PC to the ABC Host machine via ethernet cable
- Login to the ABC Host machine as ABC, password abc

Step 1




- Using the Mouse, bring up the Taskbar.
- Left click on the Folder Icon to bring up File Explorer.
- Navigate the Left-hand panel and locate your USB Device.
- Locate the saved "File Name".zip, right click and select "Extract All..."

Step 2




- Select "Browse..."
- Using the left-hand panel to locate the C: Drive.
- Navigate the "Users" folder, "Public", and select "Public Desktop"
- Then left click "Select Folder"
- Confirm the extract destination is "C:\Users\Public Desktop"

Step 3



- To test to make sure the file works.
- Double click ABC icon
- Sign in = abc
- Password = abc
- Wait for the connection to be made.

Step 4




- You are now connected to ABC Host PC.


Manually create ABC.bat file on the ABC Cleint.

Note: This is needed only if ABC.bat was not downloaded from Aktina.com.


Step 1




- At client PC, press the Windows Key and start typing "Notepad"
- Press enter or click "Open" on the right-hand side
- In Notepad type the following:



- mstsc.exe /prompt /shadow:1 /v:192.168.1.5 /admin /f /control /noconsentprompt
- note that, there is a space before the forward slash




- Once that's completed - Select "File" and then "Save As..."




- Name the File “**ABC.bat**” and save it on the desktop

Step 2



- ### Step 3





- Proceed with Installing USB over Ethernet on both the client and host PC before testing ABC Software
- Close up the remote desktop window by clicking the X on the upper right.

Installing USB over Ethernet on the Client Side

Install setup:

- You will need a license key for the USB over Ethernet install. Email support@aktina.com for a permanent product key (a 15 days trial period is in effect at installation).
- Connect the client PC to the ABC Host machine via ethernet cable
- Connect the Gating Kit device to the client PC
- Login to the ABC Host machine as ABC, password abc


Step 1





- At client PC, press the Windows Key and start typing "Edge"
- Press enter or click "Open" on the right-hand side

- While connected to the Internet or downloading from another PC with internet access. Go to the following URL link: <http://www.usb-over-ethernet.com/downloads/usb-over-ethernet.zip>
- Once typed into a web browser you should automatically download a .ZIP file
- Select "Open file" once its download


Step 2




- Once open, Double click "Setup"
- Select "I accept the agreement"
- Follow the prompts to install the software




- Continue to follow the prompts to install the software



- Continue to follow the prompts to install the software
- Then click "Next" while leaving the settings to default




- Click "Install"




- Until you get to “Finish”

Step 3




- Once the application starts.
- Locate the “Elekta Response Gating System (Communication Modem Device)”
- Please make sure the Gating kit is connected to the Client’s Laptop
- Right click Elekta Response Gating System and Select “Share”

Step 4




- Finish up by Selecting “Settings” on the top.
- Select “Preferences...”




- Under “General” tab make sure the first 4 boxes are selected and the last one turned off.
- Select “OK” once you finished.


Step 5



- With the application still open
- Click on “Help” and select “**Enter License Key...**”
- This will open a window that will allow you to copy the Product or manually type it in




- Email support@aktina.com for a permanent product key (a 15 days trial period is in effect at installation).
- Once entered you can Select “OK”




- If done correctly you will get a Prompt asking to restart the application
- We recommend doing a restart on the device and log into the account that will mainly be used


Step 6





- Bring up the application once more
- You may check to make sure its activated if logged into another profile
- Navigate to "Help" on the top and select "About"
- Close out by clicking "OK"



Step 7



- Lasty you can disable the sharing of any other device to mitigate confusion
- Right click on the device you wish not to share and select “Properties”
- Make sure you check the box for “Deny manual and automatic sharing of this device”
- You’ll see a red icon device to show its disabled





Installing USB over Ethernet Client on the Host

Install setup:


- Connect the client PC to the ABC Host machine via ethernet cable
- Connect the Gating Kit device to the client PC
- Login to the ABC Host machine as ABC Administrator, password elekta

Step 1




- At the ABC host, press the Windows Key and start typing “Edge”
- Press enter or click “Open” on the right-hand side
- While connected to the Internet or downloading from another PC with internet access. Go to the following URL link: <http://www.usb-over-ethernet.com/downloads/usb-over-ethernet-client.zip>
- Once typed into a web browser you should automatically download a .ZIP file
- Select “Open file” once its download

Step 2





- Double click on “setup_client”
- Then click on “Run”

Step 3





- Select "I accept the agreement"
- Select "Next"






- Leave the rest as default and select "Next"






- Click "Install"




- Until you get to "Finish"

Step 4




- Select "No" as this unit will not have internet access available

Step 5




- Select the PC with the Plus on the upper left

Step 6




- Type in the IP address of the ABC client device “192.168.1.55”
- Click “OK”

Step 7




- If you configured your Client Device Proceed go to the next steps. If not jump back to the following section “**Installing USB over Ethernet on the Client Side**” to install the Host software on the client side


Step 8




- Right click on “Elekta Response Gating System (Communication modem device)”
- Select “Connect”. You may need to reconnect when you sign into ABC account for the first time of software being installed. After it will automatically reconnect.
- Once you see the device connected. Right click once more and select “Properties...”




Step 9




- Under Device Connection Mode, select “Connect automatically when it is available on server”
- Select “OK”



Step 10



- On the upper left click on “Settings” and select “Preferences...”
- Check the first 4 boxes and leave the last one unchecked.
- Please do the same when logged into the ABC account.
- Once completed, you can hit “OK” and close all the application.



- On the upper left click on “Settings” and select “Preferences...”
- Check the first 4 boxes and leave the last one unchecked.
- Once completed, you can close all the application.