



Quality Control for your Active Breathing Coordinator (ABC®)

Know exactly how your ABC is performing, each and every time.

SpiroCheck Duo is a completely automated quality assurance system for the Active Breathing Coordinator. It takes the uncertainty out of your ABC performance and provides unparalleled full system calibration within two minutes. Accurate to within 1.25%, you can rest assured knowing that you are providing the optimal patient care when utilizing the ABC unit for breath holds.

Advanced firmware and a known volume of air in the automated syringe are behind the simple pass/fail test result which is readily displayed on the integrated LCD screen. Calibrated volume displacement is 1 – 3 L. PC connectivity allows users to create and print QC reports. The SpiroCheck Duo is compatible with ABC versions 2 and 3.

- Accurate to within 1.25%
- Performs full calibration within 2 minutes
- Compatible with ABC® versions 2 and 3
- PC connectivity to create and print QC reports
- Performs balloon valve and spirometry QC
- Completely automated, LCD screen displays PASS/FAIL
- Calibrated volume displacement (1 - 3 L)

